

product data



# pressure mat

## What is the Pressure Mat?

On detection of pressure exerted on the mat, an alarm output will be provided in two forms. This can be hard wired into the home equipment or wireless when used in conjunction with a Radio Output Module (ROM) or Universal Sensor.

The type of alarm call raised is dependent on how the equipment is configured.

The Pressure Mat is compatible with Lifeline home units and other Tunstall telecare enabled systems.



## Why is it needed?

The Pressure Mat can be used in the following ways:

- **Inactivity** - discreetly monitors a user's movements, eg placed under carpets and raises an alarm if no activity has been detected.
- **Activity** - monitors for activity as part of the intruder functionality available on Tunstall home units.

## Who is it for?

The Pressure Mat is for people:

- who live in a property at risk from intrusion
- who require inactivity monitoring, such as those at risk of falling

All the reassurance you need

**Tunstall**

# Features and benefits

- **Discreet installation** - with minimum visual impact
- **Dual use inactivity/intruder monitoring** - providing two types of monitoring with one device
- **Wireless operation** - the Radio Output Module can be used to provide wireless radio messages

## Technical details

Weight: 450g  
Dimensions: 700 x 400 x 10mm (WxHxD)  
Output type: normally open

### Standards:

Design, Manufacture,  
Installation and Service: IS09001: 2000

**Part Number:** *Pressure Mat GS140*  
*Radio Output Module 67005/35*

# Tunstall

[www.emergencyresponse.ie](http://www.emergencyresponse.ie)

Tunstall is a sunder member of the Continua Health Alliance

Our policy of continual development means that product specification and appearance may change without notice. Tunstall does not accept responsibility for any errors and omissions contained within this document.

© 2007 Tunstall Group Ltd. @TUNSTALL, AMIE and LIFELINE are registered trademarks.

Emergency Response Ltd Ryland Road Bunclody Enniscorthy Co Wexford Tel:  
1850 247 999 [Email: sales@emergencyresponse.ie](mailto:sales@emergencyresponse.ie)

  
Continua  
HEALTH ALLIANCE

