



Introducing the
Lifeline Vi

Vi+



Home / Away



Cancel



Help





1950 Mum, look! I've caught a jellyfish!



1960 My first date, with my future husband!



1970 I can't wait for us to be together.

Life is a journey

“ **We wanted the Lifeline Vi to reflect the time we've spent working with our customers to understand their needs for a telecare service.** ”

The Lifeline Vi is a huge leap forward to incorporate as much of their 'wish list' as we practically could in an economically viable package, and a smaller, more attractive unit. It also has unique features which make it easier to use and more flexible than its predecessors.

Richard Sharp
Group Product Development Director



1980 Let's see who can make the biggest sandcastle.



2000 A rare cuddle with my grandchildren!



TODAY There's no place like home.

The Lifeline Vi

Tunstall's sixth generation home unit has evolved. The most technically advanced, flexible and simple to use telecare platform we've ever made.

And the **best looking.**



Product Overview

Building upon our legacy of over 50 years of research, development and innovation, the Lifeline Vi sets a new benchmark in the development of telecare solutions.

Feature rich, yet easier than ever to use, and with unsurpassed attention to detail.

Clever, not complicated.



Key features

- Table stand or easy wall mount options
 - Automated regular pendant test reminder
 - Removable speaker cover for cleaning
 - Embossed cancel button for visually impaired
 - Improved speech/voice quality with automatic gain control
 - Integral ambient temperature sensor with programmable high and low alerts
- The portable installer keypad is compact and robust, and connects to the Lifeline Vi/Vi+ to enable easy programming. The keypad can also act as an ID card holder, with two attachment points allowing it to be attached to a lanyard or a belt clip. What's more, we've also made it backward compatible with the Lifeline Connect and Connect+

Enhancing your **Service**

Lifeline Vi+

The **Lifeline Vi+** is the ultimate telecare platform. As well as the suave good looks and innovative new features that make the Lifeline Vi such a winner, the Vi+ has extra functionality designed to make it even easier for you to provide an outstanding telecare service.

Improved virtual sensor management

- Not in bed/not out of bed and weekday/weekend settings
- Quick code for changing virtual bed occupancy sensor absence time

Calendar-based reminder functionality

- Single reminder, daily reminder, weekly reminder or monthly reminder
- Configurable acknowledgement of reminders – cancel button or pendant pressed to stop reminder, or no acknowledgement required option

Flexible voice recordings

- Up to 6 flexible messages with total recording time of 60 seconds (e.g. 6 x 10 seconds, 3 x 20 seconds)
- Play recorded message when a telecare sensor is activated

Feature Review

The Lifeline Vi has a wealth of features, including the following highlights ▶

Please see the Lifeline Vi solutions sheet at tunstall.com for a comprehensive feature list and full specification.

| Feature | Vi | Vi+ | |
|---|----|-----|--|
| Tunstall patented STMF technology | ✓ | ✓ | |
| Integral ambient temperature sensor | ✓ | ✓ | |
| Automated, regular pendant test reminder | ✓ | ✓ | |
| Local programming using telephone keypad | ✓ | ✓ | |
| Advanced local programming using installer keypad | ✓ | ✓ | |
| EN300 220-2 (2010) Category 1 radio receiver | ✓ | ✓ | |
| Radio interference monitor in-line with EN50134 | ✓ | ✓ | |
| BS8521 compatible | ✓ | ✓ | |
| Tunstall patented intelligent speech switching | ✓ | ✓ | |
| Plug & play telecare sensor registration | ✓ | ✓ | |
| 75% more energy efficient power adaptor | ✓ | ✓ | |
| Improved battery back up | ✓ | ✓ | |
| Telecare sensor inputs | 12 | 35 | |
| Improved virtual sensor management functionality | | ✓ | |
| Calendar-based reminder functionality | | ✓ | |
| Flexible voice recordings | | ✓ | |



sales@emergencyresponse.ie

Tunstall Emergency Response Ltd

Ryland Road, Bunclody, Enniscorthy, Co. Wexford

Tel: 1850 247999

©2012 Tunstall Group Ltd.

®Tunstall and LIFELINE are registered trademarks.

